












**MYCO Medical**  
Wound Sutures

MYCO CODE	SUTURE SIZE	NEEDLE CODE	NEEDLE DESCRIPTION	SUTURE LENGTH	NEEDLE IMAGE
<b>SILK - BLACK BRAIDED 12/BOX</b>					
SK7762	3-0	C-26	½ CIR REV 16mm	18"	X-8 
SK632	3-0	C-31	½ CIR REV 23mm	18"	X-1 
SK7772	3-0	C-6	3/8 CIR REV 19mm	18"	FS-2 
SK683	4-0	C-6	3/8 CIR REV 19mm	18"	FS-2 
SK682	5-0	C-6	3/8 CIR REV 19mm	18"	FS-2 
SK641	4-0	C-3	3/8 CIR PPT REV 13mm	18"	P-3 
SK640	5-0	C-3	3/8 CIR PPT REV 13mm	18"	P-3 
SK684	3-0	C-7	3/8 CIR REV 24mm	18"	FS-1 
SK681	6-0	PC-34	3/8 CIR PPT REV 16mm	18"	P-3 
SK685	2-0	C-8	3/8 CIR REV 26mm	18"	FS 
SK750	6-0	C-22	3/8 CIR REV 16mm	18"	FS-3 

CROSS REFERENCE	
LOOK CODE	ETHICON CODE
778B	n/a
785B	632G
784B	n/a
781B	683G
754B	682G
792B	641G
1265B 793B	640G
786B	684G
1264B	1675H
770B	685G
750B	n/a

