

Bringing Value, Reliability and Innovation to Healthcare...®

MYCO BLOOD

Collection & Infusion

VAKU-8 BLOOD COLLECTION & UNOLOK INFUSION SETS

Provide superior patient care without compromising personal safety. Both collection and infusion sets meet OSHA safety standards with needles that can be easily shielded during withdrawal utilizing single-handed activation.



Product Features

- Color coded wings for easy gauge identification
- Shielded needle minimizes risk of accidental needle stick
- OSHA-compliant shielded needle features single-handed activation
- Proprietary silconization process enhances patient comfort by reducing penetration force
- 100% latex free, constructed with kink-resistant tubing and preassembled multiple sample luer adapter

Manufactured from non-toxic, medical grade PVC on highly automated assembly machines. Needles are inspected by electronic vision camera systems to identify blunt, damaged or blocked cannula to ensure perfect needle points.



Available with and without safety shield.