



## Myco NOW offers a Fixed Winged Epidural Tuohy Point (epidural needle)

- MYCO's Huber Point epidural needle has a curved design which allows protection of the catheter from damage during catheterization procedures and its smooth insertion into the epidural space.
- The ultra thin wall of the cannula allows insertion of a catheter of a larger size. The gradual insertion at 1 cm interval allows confirmation of the depth of insertion.
- MYCO's Tuohy Needles are manufactured with fixed or detachable wings.
- Available in sterile and bulk non-sterile form.

### Superior Features

- Sharp point grinds
- Larger hubs for enhanced control
- Clear hubs for better CSF visibility

### Product Features

- Made using thin walled stainless steel for maximum flow rates
- Smooth but moderate insertion resistance makes it possible to easily identify tissue site
- Exact fit with the edge point of the Stylet with Cannula

