

Bringing Value, Reliability and Innovation to Healthcare...<sup>®</sup>

# MYCO SHARPS

## Safety

Single-handed scalpel blade removal, containment & disposal system



THE QCLICKSMART WAY



### Safety Features


- Hands positioned behind sharp
- Puncture resistant and leak proof container
- Single-handed use keeps hands out of harms way

### Easy to Use

- Right or left handed use
- 100 blade capacity
- Use with handles NO. 3,4,5,7 and 9
- Works with all manufacturers' surgical blades up to NO. 25

### Product Benefits

- Practical and reliable
- Eliminates injuries to downstream staff
- Device secured with universal, reusable mounting bracket
- When used properly the device promotes OSHA Compliance




1 CLICK EASY

Place blade flat onto landing pad - sharp edge to the **right**



2 CLICK SAFE

Insert scalpel handle until audible click



3 QCLICKSMART

Remove handle