



MYCO

MEDICAL

TheSharpChoice[®]


NMSDC Certified Diversity Supplier

RELI[®] Blunt Fill Needles

Reli[®] Blunt Fill Needles eliminate the use of a sharp when reconstituting medication, while maintaining compatibility with most syringes available in the market.

Note: Reli[®] Blunt Fill Needles are not for skin injection.

- Angled bevel reduces coring of vial stopper
- 18G ensures rapid filling and smooth flow rate suitable for viscous medications
- 1.5" needle aids penetration of different septum barriers
- Eliminates the use of a sharp, reducing the risk of a needle-stick injury
- Compatible with most syringes available in the market
- Sterile, single-use, PVC-free, not made with natural rubber latex



BRAND (COMPANY NAME)	Angled Bevel	18G	1.5"	Red Shield
Reli [®] Blunt Fill Needle (MYCO Medical)	✓	✓	✓	✓
BD [™] Blunt Fill Needle (BD)	✓	✓	✓	✓
JELCO [®] Blunt Fill Needle (Smiths Medical)	✓	✓	✓	

Reli® Blunt Fill Needles (Sterile)

MYCO CODE	GAUGE	LENGTH	UNIT OF MEASURE
BFN18G101	18G	1"	100 ea/box
BFN18G151	18G	1.5"	100 ea/box
BFN20G101	20G	1"	100 ea/box



Angled bevel

Reli® Blunt Fill Needles with Filter (Sterile)

MYCO CODE	GAUGE	LENGTH	UNIT OF MEASURE
F-BFN18G151	18G	1.5"	100 ea/box



Reli® Blunt Fill Needles
with Filter (Sterile)



Reli® Blunt Fill Needles (Sterile)

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com
www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2017 MYCO Medical, Inc. All rights reserved.
LIT-00325 Rev 2 Eif 12.14.2017