

TheSharpChoice[®]


NMSDC Certified Diversity Supplier

RELI[®] CHIBA NEEDLES


Reli[®] Chiba Needles are the value alternative to leading brands and offer superior quality in a broad range of sizes.

- Needle gauges available: 18G to 25G
- Needle lengths available: 3.5" to 10"
- Available in sterile and non-sterile packaging

- Metric markings and echogenic stylets for improved depth visibility
- Depth stopper aids accurate placement
- Large, clear hub for enhanced tactile feel
- Exact fit between stylet and cannula bevel
- Computerized grinding ensures sharp, consistent, high-quality needles




BRAND (Company)	Large Hubs	Clear Hubs	Color-Coded Stylets	Fitted, Close-Tolerance Stylet	Echogenic Stylet
Reli [®] Chiba Needles (MYCO Medical)	✓	✓	✓	✓	✓
Chiba Point (BD)			✓	✓	
Chiba Point (Kimberly Clark)		✓	✓	✓	✓



Large, clear hub with depth stopper



Metric markings and echogenic stylets



Broad range of sizes

Sterile Reli® Chiba Needles

MYCO CODE	GAUGE	LENGTH	UNIT OF MEASURE
CHE18G351	18G	3.5"	25 ea/bx; 4 bx/cs
CHE18G601	18G	6"	25 ea/bx; 4 bx/cs
CHE18G801	18G	8"	25 ea/bx; 4 bx/cs
CHE20G351	20G	3.5"	25 ea/bx; 4 bx/cs
CHE20G601	20G	6"	25 ea/bx; 4 bx/cs
CHE20G701	20G	7"	25 ea/bx; 4 bx/cs
CHE20G801	20G	8"	25 ea/bx; 4 bx/cs
CHE22G351	22G	3.5"	25 ea/bx; 4 bx/cs
CHE22G501	22G	5"	25 ea/bx; 4 bx/cs
CHE22G601	22G	6"	25 ea/bx; 4 bx/cs
CHE22G701	22G	7"	25 ea/bx; 4 bx/cs
CHE22G801	22G	8"	25 ea/bx; 4 bx/cs
CHE22G101	22G	10"	25 ea/bx; 4 bx/cs
CHE23G601	23G	6"	25 ea/bx; 4 bx/cs
CHE23G801	23G	8"	25 ea/bx; 4 bx/cs
CHE25G601	25G	6"	25 ea/bx; 4 bx/cs
CHE25G801	25G	8"	25 ea/bx; 4 bx/cs

Non-Sterile Reli® Chiba Needles

MYCO CODE	GAUGE	LENGTH	UNIT OF MEASURE
CHE18G350	18G	3.5"	100 ea/bx
CHE18G600	18G	6"	100 ea/bx
CHE18G800	18G	8"	100 ea/bx
CHE20G350	20G	3.5"	100 ea/bx
CHE20G600	20G	6"	100 ea/bx
CHE20G700	20G	7"	100 ea/bx
CHE20G800	20G	8"	100 ea/bx
CHE22G350	22G	3.5"	100 ea/bx
CHE22G500	22G	5"	100 ea/bx
CHE22G600	22G	6"	100 ea/bx
CHE22G700	22G	7"	100 ea/bx
CHE22G800	22G	8"	100 ea/bx
CHE22G100	22G	10"	100 ea/bx
CHE23G600	23G	6"	100 ea/bx
CHE23G800	23G	8"	100 ea/bx
CHE25G600	25G	6"	100 ea/bx
CHE25G800	25G	8"	100 ea/bx

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com
www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2017 MYCO Medical, Inc. All rights reserved.
LIT-00194 Rev 3 Efl 02.10.2017