



TheSharpChoice®

NMSDC Certified Diversity Supplier

DATASHEET

RELI® MULTICORE® Disposable Automatic Soft Tissue Biopsy Device

Reli® Multicore® devices are state of the art among disposable automatic biopsy instruments: they are versatile, simple to use, reliable, and feature an ergonomic design for easy cocking and firing.

- Spring-loaded stylet and cannula harvest diagnostic-quality specimens with minimal tissue displacement
- Can be operated easily with one hand, making it ideal for procedures requiring ultrasound guidance

Reli® Multicore® soft tissue biopsy devices from MYCO Medical are indicated to provide tissue samples of various soft organs and tissues, including, though not limited to, biopsies for breast, lung, liver, pancreas, kidneys, and prostate, for diagnostic sampling of abnormalities.


They are designed to provide tissue for histological examination with partial or complete removal of the imaged abnormality.

- Available in a variety of needle gauges and lengths
- Color coded for easier needle gauge identification
- Sterile, single-patient use reduces risk of infection
- Easy-to-read centimeter markings
- Echogenic tip for ease of use during imaging



Sterile Reli® Multicore® Disposable Automatic Soft Tissue Biopsy Devices

MYCO CODE	NEEDLE GAUGE	NEEDLE LENGTH	UNIT OF MEASURE
R-MLC1411	14G	11 cm	15/bx
R-MLC1415	14G	15 cm	15/bx
R-MLC1820	18G	20 cm	15/bx
R-MLC1825	18G	25 cm	15/bx



As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



MYCO
MEDICAL
TheSharpChoice®

2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com
www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2018 MYCO Medical, Inc. All rights reserved.
LIT-00526 Rev 0 Eif 11.15.2018