


TheSharpChoice[®]
NMSDC Certified Diversity Supplier

RELI[®] QUINCKE POINT Spinal Needles

Reli[®] Quincke Point Spinal Needles are the value alternative to leading brands and offer superior quality in a broad range of sizes.

- Needle gauges available: 18G to 27G
- Needle lengths available: 1.5" to 8"



- Computerized grinding ensures sharp, consistent, high-quality needles
- Exact fit between stylet and cannula bevel to minimize coring and trauma
- Large, clear hubs for enhanced tactile feel and visualization of CSF
- Smooth insertion resistance provides for easy identification of tissue layers
- Available in sterile and non-sterile packaging


BRAND (Company)	Large Hubs	Clear Hubs	Color-Coded Stylets	Fitted, Close Tolerance Stylet	Short Bevel Offerings
Reli [®] Quincke Point Spinal Needle (MYCO Medical)	✓	✓	✓	✓	✓
Spinocan (B. Braun)		✓		✓	
Quincke Spinal Needle (BD)			✓	✓	
Quincke Spinal Needle (Kimberly Clark)		✓	✓	✓	✓



Large clear hub



Matched-point grinding



Broad range of sizes

GAUGE	LENGTH	STERILE PACKAGING		NON-STERILE PACKAGING	
		MYCO Code	Unit of Measure	MYCO Code	Unit of Measure
18G	3"	SN18G301	25 ea/bx; 4 bx/cs	SN18G300	100 ea/bx
18G	3.5"	SN18G351	25 ea/bx; 4 bx/cs	SN18G350	100 ea/bx
18G	5"	SN18G501	25 ea/bx; 4 bx/cs	SN18G500	100 ea/bx
18G	6"	N/A	N/A	SN18G600	100 ea/bx
19G	3.5"	SN19G351	25 ea/bx; 4 bx/cs	SN19G350	100 ea/bx
20G	3.5"	SN20G351	25 ea/bx; 4 bx/cs	SN20G350	100 ea/bx
20G	5"	SN20G501	25 ea/bx; 4 bx/cs	SN20G500	100 ea/bx
20G	6"	SN20G601	25 ea/bx; 4 bx/cs	SN20G600	100 ea/bx
22G	1.5"	SN22G151-RW	25 ea/bx; 4 bx/cs	SN22G150-RW	100 ea/bx
22G	2.5"	SN22G251-RW	25 ea/bx; 4 bx/cs	SN22G250-RW	100 ea/bx
22G	3"	SN22G301-RW	25 ea/bx; 4 bx/cs	SN22G300-RW	100 ea/bx
22G	3.5"	SN22G351	25 ea/bx; 4 bx/cs	SN22G350	100 ea/bx
22G	4.75"	SN22G4341-RW	25 ea/bx; 4 bx/cs	SN22G4340-RW	100 ea/bx
22G	5"	SN22G501-RW	25 ea/bx; 4 bx/cs	SN22G500-RW	100 ea/bx
22G	6"	SN22G601-RW	25 ea/bx; 4 bx/cs	SN22G600-RW	100 ea/bx
22G	7"	SN22G701-RW	25 ea/bx; 4 bx/cs	SN22G700-RW	100 ea/bx
22G	8"	SN22G801-RW	25 ea/bx; 4 bx/cs	SN22G800-RW	100 ea/bx
23G	2"	N/A	N/A	SN23G200	100 ea/bx
23G	3.5"	SN23G351	25 ea/bx; 4 bx/cs	SN23G350	100 ea/bx
23G	5"	SN23G501	25 ea/bx; 4 bx/cs	SN23G500	100 ea/bx
23G	6"	SN23G601	25 ea/bx; 4 bx/cs	SN23G600	100 ea/bx
23G	7"	N/A	N/A	SN23G700	100 ea/bx
24G	3.5"	SN24G351	25 ea/bx; 4 bx/cs	SN24G350	100 ea/bx

GAUGE	LENGTH	STERILE PACKAGING		NON-STERILE PACKAGING	
		MYCO Code	Unit of Measure	MYCO Code	Unit of Measure
25G	1.5"	SN25G151	25 ea/bx; 4 bx/cs	SN25G150	100 ea/bx
25G	2.5"	SN25G251	25 ea/bx; 4 bx/cs	SN25G250	100 ea/bx
25G	3"	SN25G301	25 ea/bx; 4 bx/cs	SN25G300	100 ea/bx
25G	3.5"	SN25G351	25 ea/bx; 4 bx/cs	SN25G350	100 ea/bx
25G	5"	SN25G501	25 ea/bx; 4 bx/cs	SN25G500	100 ea/bx
25G	5"	N/A	N/A	SN25G500-RW	100 ea/bx
25G	6"	SN25G601	25 ea/bx; 4 bx/cs	SN25G600	100 ea/bx
25G	7"	SN25G701	25 ea/bx; 4 bx/cs	SN25G700	100 ea/bx
25G	8"	SN25G801	25 ea/bx; 4 bx/cs	SN25G800	100 ea/bx
26G	3.5"	SN26G351	25 ea/bx; 4 bx/cs	SN26G350	100 ea/bx
27G	3.5"	SN27G351	25 ea/bx; 4 bx/cs	SN27G350	100 ea/bx
SHORT BEVEL					
22G	3.5"	SN22G351-SB	25 ea/bx; 4 bx/cs	SN22G350-SB	100 ea/bx
22G	5"	SN22G501-SB	25 ea/bx; 4 bx/cs	SN22G500-SB	100 ea/bx
22G	6"	SN22G601-SB	25 ea/bx; 4 bx/cs	SN22G600-SB	100 ea/bx
22G	8"	SN22G801-SB	25 ea/bx; 4 bx/cs	SN22G800-SB	100 ea/bx
25G	5"	SN25G501-SB	25 ea/bx; 4 bx/cs	SN25G500-SB	100 ea/bx
METRIC MARKINGS / ECHOGENIC STYLET TIP					
18G	8"	SNME18G801	25 ea/bx; 4 bx/cs	SNME18G800	100 ea/bx
20G	3.5"	SNME20G351	25 ea/bx; 4 bx/cs	SNME20G350	100 ea/bx
22G	3.5"	SNME22G351	25 ea/bx; 4 bx/cs	SNME22G350	100 ea/bx
22G	6"	SNME22G601	25 ea/bx; 4 bx/cs	SNME22G600	100 ea/bx
25G	3.5"	SNME25G351	25 ea/bx; 4 bx/cs	SNME25G350	100 ea/bx

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to "Bringing Reliability, Innovation and Value" to the Healthcare industry with an emphasis on "physician preference" surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com
www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2017 MYCO Medical, Inc. All rights reserved.
LIT-001174 Rev 2 Eff 02.10.2017