



TheSharpChoice[®]

NMSDC Certified Diversity Supplier

UNOLOK[®] Infusion Sets

Pictured: Unolok[®]
Non-Safety Infusion Set

Unolok[®] Infusion Sets are manufactured to strict quality standards and contain unique features to improve user convenience and performance, a reliable and affordable choice for any veterinary facility

- Siliconized short bevel needle for smooth insertion
- Computerized grinding and 100% visual camera inspection ensures sharp, consistent, quality needles
- Flexible, color-coded wing for improved tactile feel and easy gauge identification
- Safety version features single-handed safety shield activation for convenience and ease of use

| BRAND (Company) | Color-Coding for Gauge Identification | Safety Option | Single-Handed Activation (Safety Version) |
|--|---------------------------------------|---------------|---|
| Unolok [®] Infusion Sets (MYCO Medical) | ✓ | ✓ | ✓ |
| Vacutainer [®] Blood Collection and Infusion Products (BD) | ✓ | ✓ | ✓ |
| Vacurette [®] Blood Collection and Infusion Products (Greiner Bio-One) | ✓ | ✓ | ✓ |
| Magellan [™] and Monoject [™] Blood Collection and Infusion Products (Kendall) | ✓ | ✓ | ✓ |
| Surshield [®] Blood Collection and Infusion Products (Terumo) | ✓ | ✓ | ✓ |

MYCO Medical's Unolok® Infusion Sets are a reliable and affordable choice for any veterinary facility, suitable for use on animals both large and small.


Sterile Unolok® Infusion Sets

| MYCO CODE | GAUGE | LENGTH | UNIT OF MEASURE |
|-----------|------------|--------|---------------------------|
| 7001-19* | 19G x 3/4" | 12" | 100/box, 10 boxes/case |
| 7001-21 | 21G x 3/4" | 12" | 100/box, 10 boxes/case |
| 7001-22* | 22G x 3/4" | 12" | 100/box, 10 boxes/case |
| 7001-23 | 23G x 3/4" | 12" | 100/box, 10 boxes/case |
| 7001-25 | 25G x 3/4" | 12" | 100/box, 10 boxes/case |


Sterile Unolok® Plus Safety Infusion Sets

| MYCO CODE | GAUGE | LENGTH | UNIT OF MEASURE |
|------------|------------|--------|---------------------------|
| 7001P-19L* | 19G x 3/4" | 8" | 100/box, 10 boxes/case |
| 7001P-21L | 21G x 3/4" | 8" | 100/box, 10 boxes/case |
| 7001P-23L | 23G x 3/4" | 8" | 100/box, 10 boxes/case |
| 7001P-25L | 25G x 3/4" | 8" | 100/box, 10 boxes/case |


**Special order with minimums and lead time*



Unolok® Infusion Sets



Color-coded wings for easy gauge identification



Unolok® Plus Safety Infusion Set

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com
www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2019 MYCO Medical, Inc. All rights reserved.
LIT-00590 Rev 0 Eff 03.20.2019