



**MYCO**  
**MEDICAL**

TheSharpChoice®


NMSDC Certified Diversity Supplier

# RELI® PRO

## ABSORBABLE SUTURES

### Package Opening Techniques

## PAPER CARD



**STEP 1:**  
Locate “TEAR”  
arrow




**STEP 2:**  
Follow “TEAR”  
path down as  
indicated




**STEP 3:**  
Open side flap to expose suture  
needle

## FOIL PACK



**STEP 3:**  
Remove LMT (less memory  
tray) from the pouch  
**NOTE:** Isopropyl-based conditioning  
fluid present

## PAPER CARD



## FOIL PACK



**STEP 4:**  
Ready to arm



**To order, learn more or receive free evaluation samples:**

Please call MYCO Medical at (800) 454-6926 or email us at [customerservice@mycomedical.com](mailto:customerservice@mycomedical.com)



2015 Production Drive  
Apex, NC 27539 USA  
[info@mycomedical.com](mailto:info@mycomedical.com)  
[www.mycomedical.com](http://www.mycomedical.com)

Some medical devices are not licensed for sale in Canada.

© 2019 MYCO Medical, Inc. All rights reserved.  
LIT-00649 Rev 0 Eff 11.25.2019